


AHSG Antibody (R35387)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R35387-100UG | 0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

Bulk quote request

| Availability | 1-3 business days |
|--------------------|--|
| Species Reactivity | Human |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 197 |
| Applications | Western Blot : 0.2-0.5ug/ml ELISA (peptide) LOD : 1:128000 |
| Limitations | This AHSG antibody is available for research use only. |



Description

Additional name(s) for this target protein: Alpha-2-HS-glycoprotein, fetuin-A

Application Notes

Optimal dilution of the AHSG antibody should be determined by the researcher.

Immunogen

Amino acids QPEGANEAVPTP were used as the immunogen for this AHSG antibody.

Storage Aliquot and store the AHSG antibody at -20oC.