

AIBZIP Antibody (Biotin Conjugate) (R35151BTN)

Catalog No.	Formulation	Size
R35151BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Biotin Conjugate
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	148327
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This AIBZIP antibody is available for research use only.



Western blot testing of human placenta lysate with biotinylated AIBZIP antibody at 1ug/ml. Predicted molecular weight: ~45kDa.

Description

Additional name(s) for this target protein: Androgen induced basic leucine zipper, CREB3L4

Application Notes

Optimal dilution of the AIBZIP antibody should be determined by the researcher.

Immunogen

Amino acids PSGRIRSVLHADEM were used as the immunogen for this AIBZIP antibody.

Storage

Aliquot and store the AIBZIP antibody at -20oC.