


AIFM1 Antibody / AIF (R33400)

Catalog No.	Formulation	Size
R33400-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9131
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This AIFM1 antibody is available for research use only.



Western blot testing of Jurkat lysate with AIFM1 antibody at 0.1ug/ml. Predicted molecular weight: ~67kDa.

Description

Additional name(s) for this target protein: Apoptosis-inducing factor, mitochondrion-associated 1

Application Notes

Optimal dilution of the AIFM1 antibody should be determined by the researcher.

1. This AIFM1 antibody detects isoforms 1, 2 & 5.

Immunogen

Amino acids DDPNVTKTLRFKQ were used as the immunogen for this AIF / AIFM1 antibody. This sequence is common to isoforms 1/2/5.

Storage

Aliquot and store the AIFM1 antibody at -20oC.