


# AIMP1 Antibody (Biotin Conjugate) (R34600BTN)

Catalog No.	Formulation	Size
R34600BTN	$0.5\ \text{mg/ml}$ in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Pig
Format	Biotin Conjugate
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q12904
Gene ID	9255
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This AIMP1 antibody is available for research use only.



### **Description**

Additional name(s) for this target protein: ARS-interacting multifunctional protein 1

## **Application Notes**

Optimal dilution of the AIMP1 antibody should be determined by the researcher.

# Immunogen

Amino acids KQQSIAGSADSKP were used as the immunogen for this AIMP1 antibody.

### **Storage**

Aliquot and store the AIMP1 antibody at -20oC.