


## AIMP1 Antibody (R34866)

Catalog No.	Formulation	Size
R34866-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q12904
Gene ID	9255
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This AIMP1 antibody is available for research use only.



Western blot testing of Jurkat lysate with AIMP1 antibody at 0.3ug/ml. Predicted molecular weight: ~34/37kDa (isoforms 1/2).

## Description

Additional name(s) for this target protein: Aminoacyl tRNA synthetase complex-interacting multifunctional protein 1, SCYE1

## Application Notes

Optimal dilution of the AIMP1 antibody should be determined by the researcher.

## **Immunogen**

Amino acids NNDVLRLEQK were used as the immunogen for this AIMP1 antibody.

## **Storage**

Aliquot and store the AIMP1 antibody at -20oC.