


## AIMP1 Antibody (R34866)

Catalog No.	Formulation	Size
R34866-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Predicted Reactivity</b>	Cow, Pig
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	Q12904
<b>Gene ID</b>	9255
<b>Applications</b>	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
<b>Limitations</b>	This AIMP1 antibody is available for research use only.



Western blot testing of Jurkat lysate with AIMP1 antibody at 0.3ug/ml. Predicted molecular weight: ~34/37kDa (isoforms 1/2).

## Description

Additional name(s) for this target protein: Aminoacyl tRNA synthetase complex-interacting multifunctional protein 1, SCYE1

## **Application Notes**

Optimal dilution of the AIMP1 antibody should be determined by the researcher.

## **Immunogen**

Amino acids NNDAVLKRLEQK were used as the immunogen for this AIMP1 antibody.

## **Storage**

Aliquot and store the AIMP1 antibody at -20oC.