


AKAP10 Antibody (R33163)

Catalog No.	Formulation	Size
R33163-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	11216
Applications	Western Blot : 0.3-1.0ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This AKAP10 antibody is available for research use only.



Western blot testing of human testis lysate with AKAP10 antibody at 0.3ug/ml. Predicted molecular weight: ~74kDa.

Description

Additional name(s) for this target protein: Mitochondrial A kinase PPKA anchor protein 10

Application Notes

Optimal dilution of the AKAP10 antibody should be determined by the researcher.

Immunogen

Amino acids AQYDQPLEKSTKL were used as the immunogen for this AKAP10 antibody.

Storage

Aliquot and store the AKAP10 antibody at -20oC.