


AKAP10 Antibody (R35661)

Catalog No.	Formulation	Size
R35661-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	11216
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This AKAP10 antibody is available for research use only.



Western blot testing of rat brain lysate with AKAP10 antibody at 1ug/ml. Predicted molecular weight: ~74kDa.

Description

Additional name(s) for this target protein: A kinase (PRKA) anchor protein 10

Application Notes

Optimal dilution of the AKAP10 antibody should be determined by the researcher.

Immunogen

Amino acids RKVKGKEQEKTS were used as the immunogen for this AKAP10 antibody.

Storage

Aliquot and store the AKAP10 antibody at -20°C.