


AKAP3 Antibody (R36240)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R36240-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|----------------------|--|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Predicted Reactivity | Cow, Dog, Mouse, Pig, Rat |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 10566 |
| Applications | Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000 |
| Limitations | This AKAP3 antibody is available for research use only. |



Western blot testing of human testis lysate with AKAP3 antibody at 1ug/ml. Predicted molecular weight: ~95kDa.

Description

Additional name(s) for this target protein: A-kinase anchor protein 3; sperm oocyte binding protein 1 (SOB1)

Application Notes

Optimal dilution of the AKAP3 antibody should be determined by the researcher.

Immunogen

Amino acids TPLQLLDWLMVNL were used as the immunogen for this AKAP3 antibody.

Storage

Aliquot and store the AKAP3 antibody at -20oC.