

AKAP6 Antibody (R34745)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R34745-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|----------------------|--|
| Availability | 1-3 business days |
| Species Reactivity | Human, Mouse |
| Predicted Reactivity | Cow, Dog, Pig, Rat |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 9472 |
| Applications | Western Blot : 1-3ug/ml (1) ELISA (peptide) LOD : 1:16000 |
| Limitations | This AKAP6 antibody is available for research use only. |



Description

Additional name(s) for this target protein: A kinase anchor protein 6, ADAP6, PRKA6, mAKAP, ADAP100

Application Notes

Optimal dilution of the AKAP6 antibody should be determined by the researcher.

1. Reported to showed a 190kDa band in lysates of mouse heart and of mouse parotid (salivary) gland.

Immunogen

Amino acids LTMSVTLSPLRSQ were used as the immunogen for this AKAP6 antibody.

Storage

Aliquot and store the AKAP6 antibody at -20oC.