

## AKAP6 Antibody (R34745)

Catalog No.	Formulation	Size
R34745-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9472
Applications	Western Blot : 1-3ug/ml (1) ELISA (peptide) LOD : 1:16000
Limitations	This AKAP6 antibody is available for research use only.



### Description

Additional name(s) for this target protein: A kinase anchor protein 6, ADAP6, PRKA6, mAKAP, ADAP100

### Application Notes

Optimal dilution of the AKAP6 antibody should be determined by the researcher.

1. Reported to showed a 190kDa band in lysates of mouse heart and of mouse parotid (salivary) gland.

## **Immunogen**

Amino acids LTMSVTLSPQRSQ were used as the immunogen for this AKAP6 antibody.

## **Storage**

Aliquot and store the AKAP6 antibody at -20oC.