


## AKT1 Antibody (F50141)

Catalog No.	Formulation	Size
F50141-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F50141-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Predicted Reactivity</b>	Mouse, Rat
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	P31749
<b>Applications</b>	Western Blot : 1:1000
<b>Limitations</b>	This AKT1 antibody is available for research use only.



AKT1 antibody western blot analysis in (1) Jurkat, (2) A375, and (3) Y79 lysate

## Description

The serine-threonine protein kinase encoded by the AKT1 gene is catalytically inactive in serum-starved primary and immortalized fibroblasts. AKT1 and the related AKT2 are activated by platelet-derived growth factor. The activation is rapid and specific, and it is abrogated by mutations in the pleckstrin homology domain of AKT1. It was shown that the activation occurs through phosphatidylinositol 3-kinase. In the developing nervous system AKT is a critical mediator of growth factor-induced neuronal survival. Survival factors can suppress apoptosis in a transcription-independent manner

by activating the serine/threonine kinase AKT1, which then phosphorylates and inactivates components of the apoptotic machinery.

For a microarray-validated AKT1 antibody supporting high-specificity detection, see our [AKT1 antibody \(clone AKT1/2552\)](#).

## Application Notes

Titration of the AKT1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A portion of amino acids 107-135 from the human protein was used as the immunogen for this AKT1 antibody.

## Storage

Aliquot the AKT1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.