


# Alcohol dehydrogenase Antibody (R36079)

Catalog No.	Formulation	Size
R36079-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	124
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : suitable ELISA (peptide) LOD : 1:8000
Limitations	This Alcohol dehydrogenase antibody is available for research use only.



#### **Description**

Additional name(s) for this target protein: Alcohol dehydrogenase 1A (ADH1A), alcohol dehydrogenase 1B (ADH1B), alcohol dehydrogenase 1C (ADH1C)

### **Application Notes**

Optimal dilution of the Alcohol dehydrogenase antibody should be determined by the researcher.

1. This Alcohol dehydrogenase antibody will detect the alpha, beta and gamma isoforms of the protein.

# **Immunogen**

Amino acids STAGKVMKCKA were used as the immunogen for this ADH1A/B/C antibody.

## **Storage**

Aliquot and store the Alcohol dehydrogenase antibody at -20oC.