


ALDH1B1 Antibody (R33498)

Catalog No.	Formulation	Size
R33498-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	219
Localization	Cytoplasmic, granular
Applications	Western Blot : 2-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This ALDH1B1 antibody is available for research use only.



Western blot testing of human liver lysate with ALDH1B1 antibody. Predicted molecular weight: ~57 kDa.

Description

Additional name(s) for this target protein: Aldehyde dehydrogenase 1 family, member B1

Application Notes

Optimal dilution of the ALDH1B1 antibody should be determined by the researcher.

Immunogen

Amino acids DKEQFERVLGYIQ were used as the immunogen for this ALDH1B1 antibody.

Storage

Aliquot and store the ALDH1B1 antibody at -20oC.