


## ALK1 Antibody (R34483)

Catalog No.	Formulation	Size
R34483-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	94
<b>Applications</b>	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
<b>Limitations</b>	This ALK1 antibody is available for research use only.



Western blot testing of MCF7 lysate with ALK1 antibody at 2ug/ml. Predicted molecular weight ~60 kDa.

### Description

Additional name(s) for this target protein: Activin receptor-like kinase 1; ACVRL1. This protein is different than the similarly named protein ALK (anaplastic lymphoma kinase), Gene ID 238.

### Application Notes

Optimal dilution of the ALK1 antibody should be determined by the researcher.

## **Immunogen**

Amino acids KISNSPEKPKVIQ were used as the immunogen for this ALK1 antibody.

## **Storage**

Aliquot and store the ALK1 antibody at -20oC.