


## ALK1 Antibody (R34483)

Catalog No.	Formulation	Size
R34483-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	94
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This ALK1 antibody is available for research use only.



Western blot testing of MCF7 lysate with ALK1 antibody at 2ug/ml. Predicted molecular weight ~60 kDa.

## Description

Additional name(s) for this target protein: Activin receptor-like kinase 1; ACVRL1. This protein is different than the similarly named protein ALK (anaplastic lymphoma kinase), Gene ID 238.

## Application Notes

Optimal dilution of the ALK1 antibody should be determined by the researcher.

## Immunogen

Amino acids KISNSPEKPKVIQ were used as the immunogen for this ALK1 antibody.

## **Storage**

Aliquot and store the ALK1 antibody at -20oC.