


AMHR2 Antibody (R34253)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R34253-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|--------------------|---|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 269 |
| Applications | Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000 |
| Limitations | This AMHR2 antibody is available for research use only. |



Western blot testing of human ovary lysate with AMHR2 antibody at 2ug/ml. Predicted molecular weight: ~63/52kDa (isoforms 1/3).

Description

Additional name(s) for this target protein: AMH type II receptor, Mullerian hormone receptor, type II

Application Notes

Optimal dilution of the AMHR2 antibody should be determined by the researcher.

1. This AMHR2 antibody is specific for isoforms 1 and 3.

Immunogen

Amino acids HPQESHFPFESCPR were used as the immunogen for this AMHR2 antibody. (1)

Storage

Aliquot and store the AMHR2 antibody at -20oC.