


ANGPTL6 Antibody (F50753)

Catalog No.	Formulation	Size
F50753-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F50753-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	Q8NI99
Applications	Western Blot : 1:1000
Limitations	This ANGPTL6 antibody is available for research use only.



Description

ANGPTL6 may play a role in the wound healing process. It may promote epidermal proliferation, remodeling and regeneration and may promote the chemotactic activity of endothelial cells and induce neovascularization. ANGPTL6 may counteract high-fat diet-induced obesity and related insulin resistance through increased energy expenditure.

Application Notes

Titration of the ANGPTL6 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 206-236 from the human protein was used as the immunogen for this ANGPTL6 antibody.

Storage

Aliquot the ANGPTL6 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.