


Anillin Antibody (R36321)

Catalog No.	Formulation	Size
R36321-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	54443
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Anillin antibody is available for research use only.



Western blot testing of mouse brain lysate with Anillin antibody at 1ug/ml. Predicted molecular weight: ~124 kDa, observed here at ~100kDa.

Description

Additional name(s) for this target protein: ANLN; Scraps

Application Notes

Optimal dilution of the Anillin antibody should be determined by the researcher.

Immunogen

Amino acids TYPDDEKRKNP were used as the immunogen for this Anillin antibody.

Storage

Aliquot and store the Anillin antibody at -20oC.