


ANKK1 Antibody (R36052)


Catalog No.	Formulation	Size
R36052-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	255239
Applications	Western Blot : 0.1-0.3ug/ml IHC (FFPE) : 5ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This ANKK1 antibody is available for research use only.



Western blot testing of human substantia nigra lysate with ANKK1 antibody at 0.1ug/ml. Predicted molecular weight: ~85 kDa. Both observed bands are blocked by addition of immunizing peptide.



IHC testing of FFPE human prostate tissue with ANKK1 antibody at 5ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Additional name(s) for this target protein: Ankyrin repeat and kinase domain containing 1

Application Notes

Optimal dilution of the ANKK1 antibody should be determined by the researcher.

Immunogen

Amino acids SRQADPNLHEAEGKT were used as the immunogen for this ANKK1 antibody.

Storage

Aliquot and store the ANKK1 antibody at -20oC.