


Annexin A11 Antibody (R36391)


Catalog No.	Formulation	Size
R36391-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	311
Applications	Western Blot : 0.2-1ug/ml Immunohistochemistry (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Annexin A11 antibody is available for research use only.



Western blot testing of HeLa lysate with Annexin A11 antibody at 0.2ug/ml. Predicted molecular weight: ~54kDa.



IHC staining of FFPE human placenta with Annexin A11 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.



Western blot testing of A) HeLa and B) Jurkat cell lysate with Annexin A11 antibody.
Predicted molecular weight: ~54kDa.

Description

Additional name(s) for this target protein: ANXA11

Application Notes

Optimal dilution of the Annexin A11 antibody should be determined by the researcher.

Immunogen

Amino acids YPGYPPPPGGYPP were used as the immunogen for this Annexin A11 antibody.

Storage

Aliquot and store the Annexin A11 antibody at -20oC.