


Anti-FOXP3 Antibody (F46753)

Catalog No.	Formulation	Size
F46753-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F46753-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Primate
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q9BZS1
Applications	Western Blot : 1:1000
Limitations	This anti-FOXP3 antibody is available for research use only.



Description

Forkhead box P3 is a member of the forkhead/winged-helix family of transcriptional regulators. Defects in this gene are the cause of immunodeficiency polyendocrinopathy, enteropathy, X-linked syndrome (IPEX), also known as X-linked autoimmunity-immunodeficiency syndrome. Alternatively spliced transcript variants encoding different isoforms have been identified.

Application Notes

Titration of the anti-FOXP3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 397-423 from the human protein was used as the immunogen for this anti-FOXP3 antibody.

Storage

Aliquot the anti-FOXP3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.

Alternate Names

Scurfin