


## AOC3 Antibody (R33411)

Catalog No.	Formulation	Size
R33411-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8639
Applications	Immunohistochemistry (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This AOC3 antibody is available for research use only.



IHC testing of FFPE human liver with AOC3 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

### Description

Additional name(s) for this target protein: Amine oxidase, copper containing 3

### Application Notes

Optimal dilution of the AOC3 antibody should be determined by the researcher.

## **Immunogen**

Amino acids DEDPSFYSA<sup>1</sup>DSI<sup>2</sup>Y were used as the immunogen for this AOC3 antibody.

## **Storage**

Aliquot and store the AOC3 antibody at -20oC.