

APOA1 Antibody (R33650)

Catalog No.	Formulation	Size
R33650-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	335
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This APOA1 antibody is available for research use only.



Western blot testing of liver lysate with APOA1 antibody at 0.5ug/ml. Predicted molecular weight ~30kDa.

Description

Additional name(s) for this target protein: Apolipoprotein A-I

Application Notes

Optimal dilution of the APOA1 antibody should be determined by the researcher.

Immunogen

Amino acids DSVTSTFSKLREQ were used as the immunogen for this APOA1 antibody.

Storage

Aliquot and store the APOA1 antibody at -20oC.