


Apoa4 Antibody (R35257)

| Catalog No. | Formulation | Size |
|--------------|--|--------|
| R35257-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$ | 100 ug |

Bulk quote request

| Availability | 1-3 business days |
|----------------------|--|
| Species Reactivity | Mouse |
| Predicted Reactivity | Rat |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 11808 |
| Applications | Immunohistochemistry (FFPE) : 2-4ug/ml ELISA (peptide) LOD : 1:32000 |
| Limitations | This Apoa4 antibody is available for research use only. |



IHC testing of FFPE mouse colon with Apoa4 antibody at 2ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Apolipoprotein A-IV

Application Notes

Optimal dilution of the Apoa4 antibody should be determined by the researcher.

Immunogen

Amino acids QETERVKEEIKKELED were used as the immunogen for this Apoa4 antibody. **Storage** Aliquot and store the Apoa4 antibody at -20oC.