


Apoa4 Antibody (R35257)

Catalog No.	Formulation	Size
R35257-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	11808
Applications	Immunohistochemistry (FFPE) : 2-4ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Apoa4 antibody is available for research use only.



IHC testing of FFPE mouse colon with Apoa4 antibody at 2ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Apolipoprotein A-IV

Application Notes

Optimal dilution of the Apoa4 antibody should be determined by the researcher.

Immunogen

Amino acids QETERVKEEIKKELED were used as the immunogen for this Apoa4 antibody.

Storage

Aliquot and store the Apoa4 antibody at -20oC.