


APOA5 Antibody (R34289)

Catalog No.	Formulation	Size
R34289-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	NP_443200.2
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This APOA5 antibody is available for research use only.



Western blot testing of human liver lysate with APOA5 antibody at 0.3ug/ml. Predicted molecular weight: ~41kDa.

Description

Additional name(s) for this target protein: Apolipoprotein A-V

Application Notes

Optimal dilution of the APOA5 antibody should be determined by the researcher.

Immunogen

Amino acids HSLHDQGHSHLGDP were used as the immunogen for this APOA5 antibody.

Storage

Aliquot and store the APOA5 antibody at -20oC.