


# APOA5 Antibody (R34289)

| Catalog No.  | Formulation   | Size   |
|--------------|---|--------|
| R34289-100UG | $0.5\ \text{mg/ml}$ in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

### **Bulk quote request**

| Availability       | 1-3 business days                                       |
|--------------------|---|
| Species Reactivity | Human   |
| Format             | Antigen affinity purified                               |
| Clonality          | Polyclonal (goat origin)                                |
| Isotype            | Goat Ig   |
| Purity             | Antigen affinity  |
| Gene ID            | NP_443200.2   |
| Applications       | Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000 |
| Limitations        | This APOA5 antibody is available for research use only. |



## **Description**

Additional name(s) for this target protein: Apolipoprotein A-V

### **Application Notes**

Optimal dilution of the APOA5 antibody should be determined by the researcher.

### **Immunogen**

Amino acids HSLHDQGHSHLGDP were used as the immunogen for this APOA5 antibody.

# **Storage** Aliquot and store the APOA5 antibody at -20oC.