


APOBEC1 Antibody (R33666)

Catalog No.	Formulation	Size
R33666-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	339
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This APOBEC1 antibody is available for research use only.



Western blot testing of mouse spleen lysate with APOBEC1 antibody at 1ug/ml.
Predicted molecular weight: ~28kDa.

Description

Additional name(s) for this target protein: Apolipoprotein B mRNA editing enzyme, catalytic polypeptide 1

Application Notes

Optimal dilution of the APOBEC1 antibody should be determined by the researcher.

Immunogen

Amino acids DVFYDPRELRKEAC were used as the immunogen for this APOBEC1 antibody.

Storage

Aliquot and store the APOBEC1 antibody at -20oC.