


## APOBEC2 Antibody (R33695)

Catalog No.	Formulation	Size
R33695-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Predicted Reactivity</b>	Mouse, Rat
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	10930
<b>Applications</b>	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
<b>Limitations</b>	This APOBEC2 antibody is available for research use only.



Western blot testing of human heart lysate with APOBEC2 antibody at 1ug/ml with [B] and without [A] blocking peptide. Predicted molecular weight: ~26kDa.

### Description

Additional name(s) for this target protein: Apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 2

### Application Notes

Optimal dilution of the APOBEC2 antibody should be determined by the researcher.

## **Immunogen**

Amino acids EAASQNGEDLENLDD were used as the immunogen for this APOBEC2 antibody.

## **Storage**

Aliquot and store the APOBEC2 antibody at -20oC.