


APOBEC3A Antibody (R36292)

Catalog No.	Formulation	Size
R36292-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	200315
Applications	Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This APOBEC3A antibody is available for research use only.



IHC testing of FFPE human tonsil with APOBEC3A antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 3A; Phorbolin 1

Application Notes

Optimal dilution of the APOBEC3A antibody should be determined by the researcher.

Immunogen

Amino acids TSNFNNGIGRHKTY were used as the immunogen for this APOBEC3A antibody.

Storage

Aliquot and store the APOBEC3A antibody at -20oC.