


## APOBEC3A Antibody (R36292)

Catalog No.	Formulation	Size
R36292-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	200315
Applications	Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This APOBEC3A antibody is available for research use only.



IHC testing of FFPE human tonsil with APOBEC3A antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

### Description

Additional name(s) for this target protein: Apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 3A; Phorbolin 1

### Application Notes

Optimal dilution of the APOBEC3A antibody should be determined by the researcher.

### Immunogen

Amino acids TSNFNNGIGRHKTY were used as the immunogen for this APOBEC3A antibody.

## **Storage**

Aliquot and store the APOBEC3A antibody at -20oC.