


## APOBEC3C Antibody (R35672)

Catalog No.	Formulation	Size
R35672-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	27350
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This APOBEC3C antibody is available for research use only.



Western blot testing of Daudi lysate with APOBEC3C antibody at 0.5ug/ml. Predicted molecular weight: ~23kDa.

## Description

Additional name(s) for this target protein: Apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 3C

## Application Notes

Optimal dilution of the APOBEC3C antibody should be determined by the researcher.

## Immunogen

Amino acids RLLKRRLRESLQ were used as the immunogen for this APOBEC3C antibody.

## **Storage**

Aliquot and store the APOBEC3C antibody at -20oC.