


APOBEC3D Antibody (R35050)

Catalog No.	Formulation	Size
R35050-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	140564
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This APOBEC3D antibody is available for research use only.



Western blot testing of A549 lysate with APOBEC3D antibody at 0.3ug/ml. Predicted molecular weight: ~46kDa.

Description

Additional name(s) for this target protein: Apolipoprotein B mRNA editing enzyme catalytic polypeptide-like 3D

Application Notes

Optimal dilution of the APOBEC3D antibody should be determined by the researcher.

Immunogen

Amino acids PKRQSNHRQEVYFR were used as the immunogen for this APOBEC3D antibody.

Storage

Aliquot and store the APOBEC3D antibody at -20oC.