


APOBEC3G Antibody (R35950)

Catalog No.	Formulation	Size
R35950-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	60489
Localization	Cytoplasmic, membranous
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This APOBEC3G antibody is available for research use only.



IHC testing of FFPE human kidney with APOBEC3G antibody at 4ug/ml. HIER: steamed with pH9 Tris/EDTA buffer, HRP-staining.

Description

Additional name(s) for this target protein: Apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 3G; ARP9

Application Notes

Optimal dilution of the APOBEC3G antibody should be determined by the researcher.

Immunogen

Amino acids SKWRKLHRDQE were used as the immunogen for this APOBEC3G antibody.

Storage

Aliquot and store the APOBEC3G antibody at -20oC.