


APOC1 Antibody (R34014)

Catalog No.	Formulation	Size
R34014-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	341
Applications	Immunohistochemistry (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This APOC1 antibody is available for research use only.



IHC testing of FFPE human liver with APOC1 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Apolipoprotein C-I

Application Notes

Optimal dilution of the APOC1 antibody should be determined by the researcher.

Immunogen

Amino acids ETFQKVKEKLKIDS were used as the immunogen for this APOC1 antibody.

Storage

Aliquot and store the APOC1 antibody at -20oC.