


APOH Antibody / Apolipoprotein H (R35271)

Catalog No.	Formulation	Size
R35271-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	350
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This APOH antibody is available for research use only.



Western blot testing of human albumin-depleted plasma lysate with APOH antibody at 1ug/ml. Predicted molecular weight: ~38/50 kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: Apolipoprotein H

Application Notes

Optimal dilution of the APOH antibody should be determined by the researcher.

Immunogen

Amino acids QGERVKIQEKFKN were used as the immunogen for this APOH antibody.

Storage

Aliquot and store the APOH antibody at -20oC.