


APOL2 Antibody (R36192)

Catalog No.	Formulation	Size
R36192-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23780
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This APOL2 antibody is available for research use only.



Western blot testing of human frontal cortex lysate with APOL2 antibody at 0.03ug/ml.
Predicted molecular weight: ~38kDa.

Description

Additional name(s) for this target protein: Apolipoprotein L 2; APOL-II

Application Notes

Optimal dilution of the APOL2 antibody should be determined by the researcher.

Immunogen

Amino acids TKIHEMLQPGQD were used as the immunogen for this APOL2 antibody.

Storage

Aliquot and store the APOL2 antibody at -20oC.