


APOL3 Antibody / Apolipoprotein L3 (F53959)

Catalog No.	Formulation	Size
F53959-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F53959-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	O95236
Applications	Western Blot : 1:500-2000
Limitations	This APOL3 antibody is available for research use only.



Western blot testing of human 1) A549 and 2) HepG2 cell lysate with APOL3 antibody at 1:2000. Predicted molecular weight: 44 kDa.

Description

Apolipoprotein L3 may affect the movement of lipids in the cytoplasm or allow the binding of lipids to organelles.

Application Notes

Titration of the APOL3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 57-90 from the human protein was used as the immunogen for the APOL3 antibody.

Storage

Aliquot the APOL3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.