


APOL3 Antibody (R33879)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R33879-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|--------------------|---|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 80833 |
| Applications | Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000 |
| Limitations | This APOL3 antibody is available for research use only. |



Western blot testing of Daudi lysate with APOL3 antibody at 1ug/ml. Predicted molecular weight: ~44/37/22kDa (isoforms 1/2/3).

Description

Additional name(s) for this target protein: Apolipoprotein L3

Application Notes

Optimal dilution of the APOL3 antibody should be determined by the researcher.

1. This APOL3 antibody is specific for isoforms 1 and 2/3 of the protein.

Immunogen

Amino acids ELTQIYQRLNPCHTH were used as the immunogen for this APOL3 antibody.

Storage

Aliquot and store the APOL3 antibody at -20oC.