


APOL5 Antibody (R34959)

Catalog No.	Formulation	Size
R34959-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	80831
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This APOL5 antibody is available for research use only.



Western blot testing of human frontal cortex lysate with APOL5 antibody at 1 ug/ml.
Predicted molecular weight: ~47kDa.

Description

Additional name(s) for this target protein: Apolipoprotein L5

Application Notes

Optimal dilution of the APOL5 antibody should be determined by the researcher.

Immunogen

Amino acids PARKGRQAPGRHRQ were used as the immunogen for this APOL5 antibody.

Storage

Aliquot and store the APOL5 antibody at -20°C.