


# **APXL Antibody / SHROOM2 (R34365)**

Catalog No.	Formulation	Size
R34365-100UG	$0.5\ \text{mg/ml}$ in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	357
Applications	Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This APXL antibody is available for research use only.



IHC testing of FFPE human kidney with APXL antibody at 4ug/ml. HIER: steamed with pH9 Tris/EDTA buffer, HRP-staining.

### **Description**

Additional name(s) for this target protein: HSAPXL, Apical protein-like; SHROOM2

## **Application Notes**

Optimal dilution of the APXL antibody should be determined by the researcher.

### **Immunogen**

Amino acids KCLLDSLQPERGK were used as the immunogen for this APXL antibody. **Storage** Aliquot and store the APXL antibody at -20oC.