

ARF1 Antibody / ADP-ribosylation factor 1 (R33779)

Catalog No.	Formulation	Size
R33779-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	375; 376, 377, 378
Applications	Western Blot : 1-3ug/ml (1) ELISA (peptide) LOD : 1:8000
Limitations	This ARF1 antibody is available for research use only.



Western blot testing of HepG2 cell lysate with ARF1 antibody at 1ug/ml. Predicted molecular weight ~21kDa.

Description

Additional name(s) for this target protein: ADP-ribosylation factor 1

Application Notes

Optimal dilution of the ARF1 antibody should be determined by the researcher.

1. This Ab will detect ADP-ribosylation factor isoforms 1-4.

Immunogen

Amino acids EGLDWLSNQLRNQK were used as the immunogen for this ARF1 antibody.

Storage

Aliquot and store the ARF1 antibody at -20oC.