


ARF4 Antibody (R35883)


Catalog No.	Formulation	Size
R35883-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	378
Localization	Cytoplasmic, membranous
Applications	Western Blot : 0.3-1ug/ml IHC (FFPE) : 5ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This ARF4 antibody is available for research use only.



Western blot testing of mouse NIH3T3 [A] and human HeLa [B] lysates with ARF4 antibody at 0.5ug/ml. Predicted molecular weight: ~20 kDa.



IHC testing of FFPE human thymus tissue with ARF4 antibody at 5ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Additional name(s) for this target protein: ADP-ribosylation factor 4

Application Notes

Optimal dilution of the ARF4 antibody should be determined by the researcher.

Immunogen

Amino acids SEMTDKLGLQSLRN were used as the immunogen for this ARF4 antibody.

Storage

Aliquot and store the ARF4 antibody at -20oC.