


ARFGAP3 Antibody (R34812)

Catalog No.	Formulation	Size
R34812-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q9NP61
Gene ID	26286
Applications	Western Blot : 2-4ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This ARFGAP3 antibody is available for research use only.



Western blot testing of mouse brain with ARFGAP3 antibody at 2ug/ml. Predicted molecular weight: ~57/52kDa (isoforms 1/2).

Description

Additional name(s) for this target protein: ADP-ribosylation factor GTPase activating protein 3

Application Notes

Optimal dilution of the ARFGAP3 antibody should be determined by the researcher.

Immunogen

Amino acids NGVVTSIQDRYGS were used as the immunogen for this ARFGAP3 antibody. This sequence is common to isoforms 1 and 2 (UniProt # Q9NP61).

Storage

Aliquot and store the ARFGAP3 antibody at -20oC.