


# ARFRP1 Antibody / ADP-ribosylation factor-related protein 1 [clone 28A57] (RQ8419)

Catalog No.	Formulation	Size
RQ8419	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

#### Recombinant RABBIT MONOCLONAL

## **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	28A57
Purity	Affinity chromatography
UniProt	Q13795
Applications	Western Blot : 1:500-1:2000
Limitations	This ARFRP1 antibody is available for research use only.



Western blot testing of human HEK293 cell lysate with ARFRP1 antibody. Predicted molecular weight ~23 kDa.

#### **Description**

Trans-Golgi-associated GTPase that regulates protein sorting. [UniProt]

### **Application Notes**

Optimal dilution of the ARFRP1 antibody should be determined by the researcher.

## **Immunogen**

A synthetic peptide derived from human ARFRP1 / ARP was used as the immunogen for the ARFRP1 antibody. **Storage** Store the ARFRP1 antibody at -20oC.