


ARH Antibody (R33381)


Catalog No.	Formulation	Size
R33381-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	26119
Applications	IF/ICC : 10ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This ARH antibody is available for research use only.



IF/ICC testing of fixed and permeabilized human A431 cells with ARH antibody (green) at 10ug/ml and DAPI nuclear stain (blue).



IF/ICC testing of fixed and permeabilized human U-2 OS cells with ARH antibody (green) at 10ug/ml and DAPI nuclear stain (blue).

Description

Additional name(s) for this target protein: Autosomal recessive hypercholesterolemia protein, LDL receptor adaptor protein; LDLRAP1

Application Notes

Optimal dilution of the ARH antibody should be determined by the researcher.

Immunogen

Amino acids DALKSAGRALIRS were used as the immunogen for this ARH antibody.

Storage

Aliquot and store the ARH antibody at -20oC.