


ARHGDIG Antibody (R33721)

Catalog No.	Formulation	Size
R33721-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	398
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This ARHGDIG antibody is available for research use only.



Western blot testing of HeLa lysate with ARHGDIG antibody at 0.5ug/ml. Predicted molecular weight: ~25kDa.

Description

Additional name(s) for this target protein: Rho GDP dissociation inhibitor (GDI) gamma; RHOGDI-3; RhoGDI gamma

Application Notes

Optimal dilution of the ARHGDIG antibody should be determined by the researcher.

Immunogen

Amino acids EAVPEYRAPGRK were used as the immunogen for this ARHGDIG antibody.

Storage

Aliquot and store the ARHGDIG antibody at -20oC.