


ARHGEF4 Antibody (R35175)


Catalog No.	Formulation	Size
R35175-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	50649
Applications	Western Blot : 2-4ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This ARHGEF4 antibody is available for research use only.



Western blot testing of human brain lysate with ARHGEF4 antibody at 2ug/ml. Predicted molecular weight: ~79/76kDa (isoforms 1/2).



Western blot of HEK293 lysate overexpressing ARHGEF4 probed with ARHGEF4 antibody (mock transfection in lane 1). Predicted molecular weight: ~79/76kDa (isoforms 1/2).

Description

Additional name(s) for this target protein: APC-stimulated guanine nucleotide exchange factor 4

Application Notes

Optimal dilution of the ARHGEF4 antibody should be determined by the researcher.

Immunogen

Amino acids PWEEPAGEKPSCS were used as the immunogen for this ARHGEF4 antibody.

Storage

Aliquot and store the ARHGEF4 antibody at -20oC.