


ARIH2 Antibody (R33946)

Catalog No.	Formulation	Size
R33946-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10425
Applications	Western Blot : 1-3ug/ml Immunofluorescence : suitable Immunoprecipitation : suitable ELISA (peptide) LOD : 1:8000
Limitations	This ARIH2 antibody is available for research use only.



Western blot testing of HeLa [A] and K562 [B] nuclear lysates with ARIH2 antibody at 2ug/ml. Predicted molecular weight: ~58kDa.

Description

Additional name(s) for this target protein: Ariadne homolog 2, TRIAD1

Application Notes

Optimal dilution of the ARIH2 antibody should be determined by the researcher.

Immunogen

Amino acids EQRRRTLLKDFHDT were used as the immunogen for this ARIH2 antibody.

Storage

Aliquot and store the ARIH2 antibody at -20oC.